

Export images

1. Go to: File>Save/Export to Tiff
2. In window that opens:
 - a. Use “Folder selection” to navigate to where your images are
 - b. In the “all images” box on the left highlight the stitched images you want to export (all should end in “_001.nd2”)
 - c. Use the “Feed” button to move them to the right box.
 - a. The box on the right side is images that will be batch processed
 - d. Under “Output Folder” Check the box that says “Same as input”
 - e. In the bottom part of the window select the “OME-TIFF” option
 - f. Click “Export” to process images!

